

## Cyclic AMP-dependent transcription factor ATF-2

Catalog No: #AP80983



Package Size: #AP80983-1 50ug #AP80983-2 100ug #AP80983-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cyclic AMP-dependent transcription factor ATF-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	cAMP-dependent transcription factor ATF-2, Activating transcription factor 2, MXBP protein, cAMP response element-binding protein CRE-BP1
Accession No.	Uniprot ID:P16951
Uniprot	P16951
GeneID	11909;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Atf2

Note: This product is for in vitro research use only