

B-cell receptor-associated protein 31

Catalog No: #AP80986



Package Size: #AP80986-1 50ug #AP80986-2 100ug #AP80986-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell receptor-associated protein 31
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BCR-associated protein 31,p28,Bap31
Accession No.	Uniprot ID:Q61335
Uniprot	Q61335
GeneID	27061;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bcap31

Note: This product is for in vitro research use only