

D-beta-hydroxybutyrate dehydrogenase, mitochondrial

Catalog No: #AP80989



Package Size: #AP80989-1 50ug #AP80989-2 100ug #AP80989-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3-hydroxybutyrate dehydrogenase,BDH,Short chain dehydrogenase/reductase family 9C member 1,BDH,SDR9C1
Accession No.	Uniprot ID:Q02338
Uniprot	Q02338
GeneID	622;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BDH1

Note: This product is for in vitro research use only