

Baculoviral IAP repeat-containing protein 1

Catalog No: #AP80998



Package Size: #AP80998-1 50ug #AP80998-2 100ug #AP80998-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neuronal apoptosis inhibitory protein,BIRC1
Accession No.	Uniprot ID:Q13075
Uniprot	Q13075
GeneID	4671;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NAIP

Note: This product is for in vitro research use only