

B-cell lymphoma 3 protein homolog

Catalog No: #AP81001



Package Size: #AP81001-1 50ug #AP81001-2 100ug #AP81001-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell lymphoma 3 protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BCL-3
Accession No.	Uniprot ID:Q9Z2F6
Uniprot	Q9Z2F6
GeneID	12051;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Bcl3

Note: This product is for in vitro research use only