

Cyclic AMP-dependent transcription factor ATF-4

Catalog No: #AP81002



Package Size: #AP81002-1 50ug #AP81002-2 100ug #AP81002-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	cAMP-dependent transcription factor ATF-4,Activating transcription factor 4,C/EBP-related ATF,C/ATF,Tax-responsive enhancer element-binding protein 67 homolog,TaxREB67 homolog
Accession No.	Uniprot ID:Q06507
Uniprot	Q06507
GeneID	11911;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atf4

Note: This product is for in vitro research use only