

# Ubiquitin-like protein ATG12

Catalog No: #AP81005



Package Size: #AP81005-1 50ug #AP81005-2 100ug #AP81005-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubiquitin-like protein ATG12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Autophagy-related protein 12,APG12-like,APG12,APG12L
Accession No.	Uniprot ID:O94817
Uniprot	O94817
GeneID	9140;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ATG12

Note: This product is for in vitro research use only