

Protein atonal homolog 1

Catalog No: #AP81006



Package Size: #AP81006-1 50ug #AP81006-2 100ug #AP81006-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein atonal homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Helix-loop-helix protein mATH-1,mATH1,Ath1
Accession No.	Uniprot ID:P48985
Uniprot	P48985
GeneID	11921;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atoh1

Note: This product is for in vitro research use only