

Brain-specific serine protease 4

Catalog No: #AP81010



Package Size: #AP81010-1 50ug #AP81010-2 100ug #AP81010-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Brain-specific serine protease 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BSSP-4, Serine protease 22, Serine protease 26, Tryptase epsilon, Bssp4, Prss26
Accession No.	Uniprot ID:Q9ER10
Uniprot	Q9ER10
GeneID	70835;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prss22

Note: This product is for in vitro research use only