

Acetylserotonin O-methyltransferase

Catalog No: #AP81011



Package Size: #AP81011-1 50ug #AP81011-2 100ug #AP81011-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetylserotonin O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hydroxyindole O-methyltransferase,Hiomt
Accession No.	Uniprot ID:D3KU66
Uniprot	D3KU66
GeneID	107626;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Asmt

Note: This product is for in vitro research use only