

A disintegrin and metalloproteinase with thrombospondin motifs 1

Catalog No: #AP81013

Package Size: #AP81013-1 50ug #AP81013-2 100ug #AP81013-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ADAM-TS 1
Accession No.	Uniprot ID:Q9WUQ1
Uniprot	Q9WUQ1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adamts1

Note: This product is for in vitro research use only