

Beta-1,4-galactosyltransferase 1

Catalog No: #AP81015



Package Size: #AP81015-1 50ug #AP81015-2 100ug #AP81015-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1,4-galactosyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Beta-1,4-GalTase 1,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 1,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 1,GALT,GGTB2
Accession No.	Uniprot ID:P08037
Uniprot	P08037
GeneID	281781;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:B4GALT1

Note: This product is for in vitro research use only