

Basigin

Catalog No: #AP81016



Package Size: #AP81016-1 50ug #AP81016-2 100ug #AP81016-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Basigin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | 5A11 antigen,Blood-brain barrier HT7 antigen,Neurothelin |
| Accession No. | Uniprot ID:P17790 |
| Uniprot | P17790 |
| GeneID | 770363; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BSG

Note: This product is for in vitro research use only