

## Beta-secretase 2

Catalog No: #AP81017



Package Size: #AP81017-1 50ug #AP81017-2 100ug #AP81017-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-secretase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-site amyloid precursor protein cleaving enzyme 2, Beta-site APP cleaving enzyme 2, Memapsin-1, Membrane-associated aspartic protease 1, Theta-secretase
Accession No.	Uniprot ID:Q6IE75
Uniprot	Q6IE75
GeneID	288227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Bace2

Note: This product is for in vitro research use only