

Betaine--homocysteine S-methyltransferase 1

Catalog No: #AP81024



Package Size: #AP81024-1 50ug #AP81024-2 100ug #AP81024-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Betaine--homocysteine S-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:O35490
Uniprot	O35490
GeneID	12116;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bhmt

Note: This product is for in vitro research use only