

Carbonic anhydrase 12

Catalog No: #AP81026



Package Size: #AP81026-1 50ug #AP81026-2 100ug #AP81026-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carbonate dehydratase XII, Carbonic anhydrase XII, CA-XII, Car12
Accession No.	Uniprot ID:Q8CI85
Uniprot	Q8CI85
GeneID	76459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ca12

Note: This product is for in vitro research use only