

Betaine--homocysteine S-methyltransferase 1

Catalog No: #AP81031



Package Size: #AP81031-1 50ug #AP81031-2 100ug #AP81031-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Betaine--homocysteine S-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q5I597
Uniprot	Q5I597
GeneID	497025;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BHMT

Note: This product is for in vitro research use only