

S-methylmethionine--homocysteine S-methyltransferase BHMT2



Catalog No: #AP81032

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81032-1 50ug #AP81032-2 100ug #AP81032-3 1mg

Description

Product Name	S-methylmethionine--homocysteine S-methyltransferase BHMT2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SMM-hcy methyltransferase,Betaine--homocysteine S-methyltransferase 2
Accession No.	Uniprot ID:Q9H2M3
Uniprot	Q9H2M3
GeneID	23743;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BHMT2

Note: This product is for in vitro research use only