

Bone morphogenetic protein 7

Catalog No: #AP81034



Package Size: #AP81034-1 50ug #AP81034-2 100ug #AP81034-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMP-7,Osteogenic protein 1,OP-1,Bmp-7,Op1
Accession No.	Uniprot ID:P23359
Uniprot	P23359
GeneID	12162;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bmp7

Note: This product is for in vitro research use only