

Beta-secretase 2

Catalog No: #AP81036



Package Size: #AP81036-1 50ug #AP81036-2 100ug #AP81036-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-secretase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aspartyl protease 1,ASP1,Beta-site amyloid precursor protein cleaving enzyme 2,Beta-site APP cleaving enzyme 2,Memapsin-1,Membrane-associated aspartic protease 1,Theta-secretase
Accession No.	Uniprot ID:Q9JL18
Uniprot	Q9JL18
GeneID	56175;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bace2

Note: This product is for in vitro research use only