

Beta-1,3-galactosyltransferase 2

Catalog No: #AP81039



Package Size: #AP81039-1 50ug #AP81039-2 100ug #AP81039-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1,3-galactosyltransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta-1,3-GalTase 2,UDP-Gal:betaGlcNAc beta 1,3-galactosyltransferase-II
Accession No.	Uniprot ID:O54905
Uniprot	O54905
GeneID	26878;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:B3galt2

Note: This product is for in vitro research use only