

## Gamma Tubulin Mouse Monoclonal Antibody(7C1)

Catalog No: #44166

Package Size: #44166-1 50ul #44166-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Gamma Tubulin Mouse Monoclonal Antibody(7C1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7C1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Gamma Tubulin Mouse Monoclonal Antibody(7C1) detects endogenous levels of TUBG1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Gamma Tubulin
Other Names	TUBG1
Accession No.	Swiss-Prot:P23258Gene ID:7283
Uniprot	P23258
GeneID	7283;
Calculated MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using gamma Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using gamma Tutulin Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only