

B-cell receptor-associated protein 31

Catalog No: #AP81040



Package Size: #AP81040-1 50ug #AP81040-2 100ug #AP81040-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	B-cell receptor-associated protein 31
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BCR-associated protein 31,6C6-AG tumor-associated antigen,Protein CDM,p28,BAP31,DXS1357E
Accession No.	Uniprot ID:P51572
Uniprot	P51572
GeneID	10134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCAP31

Note: This product is for in vitro research use only