

## ATP synthase-coupling factor 6, mitochondrial

Catalog No: #AP81041

Package Size: #AP81041-1 50ug #AP81041-2 100ug #AP81041-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ATP synthase-coupling factor 6, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATPase subunit F6
Accession No.	Uniprot ID:P21571
Uniprot	P21571
GeneID	94271;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Atp5j

Note: This product is for in vitro research use only