

ATP synthase subunit beta, mitochondrial

Catalog No: #AP81042



Package Size: #AP81042-1 50ug #AP81042-2 100ug #AP81042-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP synthase subunit beta, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_6118
Accession No.	Uniprot ID:Q5ZLC5
Uniprot	Q5ZLC5
GeneID	426673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATP5B

Note: This product is for in vitro research use only