

# A disintegrin and metalloproteinase with thrombospondin motifs 13

Catalog No: #AP81043

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP81043-1 50ug #AP81043-2 100ug #AP81043-3 1mg

## Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 13,von Willebrand factor-cleaving protease,vWF-CP,C9orf8,UNQ6102/PRO20085
Accession No.	Uniprot ID:Q76LX8
Uniprot	Q76LX8
GeneID	11093;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADAMTS13

Note: This product is for in vitro research use only