

Betaine--homocysteine S-methyltransferase 1

Catalog No: #AP81045



Package Size: #AP81045-1 50ug #AP81045-2 100ug #AP81045-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Betaine--homocysteine S-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q93088
Uniprot	Q93088
GeneID	635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BHMT

Note: This product is for in vitro research use only