

Baculoviral IAP repeat-containing protein 6

Catalog No: #AP81046



Package Size: #AP81046-1 50ug #AP81046-2 100ug #AP81046-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BIR repeat-containing ubiquitin-conjugating enzyme, BRUCE, RING-type E3 ubiquitin transferase BIRC6, Ubiquitin-conjugating BIR domain enzyme apollon, APOLLON, KIAA1289
Accession No.	Uniprot ID: Q9NR09
Uniprot	Q9NR09
GeneID	57448;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: BIRC6

Note: This product is for in vitro research use only