

Asporin

Catalog No: #AP81048



Package Size: #AP81048-1 50ug #AP81048-2 100ug #AP81048-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Asporin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Periodontal ligament-associated protein 1, PLAP-1, PLAP1, SLRR1C, UNQ215/PRO241
Accession No.	Uniprot ID:Q9BXN1
Uniprot	Q9BXN1
GeneID	54829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ASPN

Note: This product is for in vitro research use only