

Serine/threonine-protein kinase ATR

Catalog No: #AP81052



Package Size: #AP81052-1 50ug #AP81052-2 100ug #AP81052-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase ATR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ataxia telangiectasia and Rad3-related protein,FRAP-related protein 1,FRP1
Accession No.	Uniprot ID:Q13535
Uniprot	Q13535
GeneID	545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATR

Note: This product is for in vitro research use only