

Protein BTG3

Catalog No: #AP81053



Package Size: #AP81053-1 50ug #AP81053-2 100ug #AP81053-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein BTG3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Abundant in neuroepithelium area protein,BTG family member 3,Protein Tob5,ANA,TOB5
Accession No.	Uniprot ID:Q14201
Uniprot	Q14201
GeneID	10950;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BTG3

Note: This product is for in vitro research use only