

A disintegrin and metalloproteinase with thrombospondin motifs 2

Catalog No: #AP81056

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81056-1 50ug #AP81056-2 100ug #AP81056-3 1mg

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM-TS 2,Procollagen I N-proteinase,PC I-NP,Procollagen I/II amino propeptide-processing enzyme,Procollagen N-endopeptidase,pNPI
Accession No.	Uniprot ID:Q8C9W3
Uniprot	Q8C9W3
GeneID	216725;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adamts2

Note: This product is for in vitro research use only