

Ataxin-10

Catalog No: #AP81057



Package Size: #AP81057-1 50ug #AP81057-2 100ug #AP81057-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ataxin-10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain protein E46 homolog, Spinocerebellar ataxia type 10 protein, SCA10
Accession No.	Uniprot ID:Q9UBB4
Uniprot	Q9UBB4
GeneID	25814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATXN10

Note: This product is for in vitro research use only