

# Phospho-Stat3 (Y705) Mouse Monoclonal Antibody(1E9)

Catalog No: #12835

Package Size: #12835-1 50ul #12835-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Phospho-Stat3 (Y705) Mouse Monoclonal Antibody(1E9)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1.00E+09
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Phospho-Stat3 (Y705) Mouse Monoclonal Antibody(1E9) detects endogenous levels of STAT3
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Stat3Mouse
Modification	Phospho
Other Names	STAT3
Accession No.	Swiss-Prot:P40763Gene ID:6774
Uniprot	P40763
GeneID	6774;
Calculated MW	79,86kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Muscle Actin Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only