

Calcitonin gene-related peptide

Catalog No: #AP81063



Package Size: #AP81063-1 50ug #AP81063-2 100ug #AP81063-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CGRP,CALC
Accession No.	Uniprot ID:P30880
Uniprot	P30880
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CALCA

Note: This product is for in vitro research use only