

Bactericidal permeability-increasing protein

Catalog No: #AP81064



Package Size: #AP81064-1 50ug #AP81064-2 100ug #AP81064-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bactericidal permeability-increasing protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BPI
Accession No.	Uniprot ID:P17453
Uniprot	P17453
GeneID	280734;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BPI

Note: This product is for in vitro research use only