

Scavenger receptor cysteine-rich type 1 protein M130

Catalog No: #AP81065



Package Size: #AP81065-1 50ug #AP81065-2 100ug #AP81065-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scavenger receptor cysteine-rich type 1 protein M130
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD163,M130
Accession No.	Uniprot ID:P85521
Uniprot	P85521
GeneID	533844;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD163

Note: This product is for in vitro research use only