

Breast cancer type 2 susceptibility protein homolog

Catalog No: #AP81066



Package Size: #AP81066-1 50ug #AP81066-2 100ug #AP81066-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Breast cancer type 2 susceptibility protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fanconi anemia group D1 protein homolog,Fancd1
Accession No.	Uniprot ID:P97929
Uniprot	P97929
GeneID	12190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Brca2

Note: This product is for in vitro research use only