

Breast cancer metastasis-suppressor 1 homolog

Catalog No: #AP81067



Package Size: #AP81067-1 50ug #AP81067-2 100ug #AP81067-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Breast cancer metastasis-suppressor 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q5M7T3
Uniprot	Q5M7T3
GeneID	108348098;293668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Brms1

Note: This product is for in vitro research use only