

B-cell lymphoma/leukemia 11A

Catalog No: #AP81069



Package Size: #AP81069-1 50ug #AP81069-2 100ug #AP81069-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell lymphoma/leukemia 11A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BCL-11A,B-cell CLL/lymphoma 11A,COUP-TF-interacting protein 1,Ecotropic viral integration site 9 protein homolog,EVI-9,Zinc finger protein 856,CTIP1,EVI9,KIAA1809,ZNF856
Accession No.	Uniprot ID:Q9H165
Uniprot	Q9H165
GeneID	53335;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCL11A

Note: This product is for in vitro research use only