

LC3B Mouse Monoclonal Antibody(9H5)

Catalog No: #44169

Package Size: #44169-1 50ul #44169-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LC3B Mouse Monoclonal Antibody(9H5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9H5
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	LC3B Mouse Monoclonal Antibody(9H5) detects endogenous levels of MAP1LC3B
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	LC3B
Other Names	MAP1LC3B
Accession No.	Swiss-Prot:Q9GZQ8Gene ID:81631
Uniprot	Q9GZQ8
GeneID	81631;
Calculated MW	14,16kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using LC3B Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only