

Bone morphogenetic protein 4

Catalog No: #AP81071



Package Size: #AP81071-1 50ug #AP81071-2 100ug #AP81071-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BMP-4,Bone morphogenetic protein 2B,BMP-2B,BMP2B,DVR4
Accession No.	Uniprot ID:P12644
Uniprot	P12644
GeneID	652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BMP4

Note: This product is for in vitro research use only