

Monofunctional C1-tetrahydrofolate synthase, mitochondrial

Catalog No: #AP81073

Package Size: #AP81073-1 50ug #AP81073-2 100ug #AP81073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monofunctional C1-tetrahydrofolate synthase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Formyltetrahydrofolate synthetase,Fthfsc1
Accession No.	Uniprot ID:Q3V3R1
Uniprot	Q3V3R1
GeneID	270685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mthfd11

Note: This product is for in vitro research use only