

Protein BTG3

Catalog No: #AP81078



Package Size: #AP81078-1 50ug #AP81078-2 100ug #AP81078-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein BTG3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rBTG3,BTG family member 3
Accession No.	Uniprot ID:O88677
Uniprot	O88677
GeneID	54230;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Btn3

Note: This product is for in vitro research use only