

Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7)

Catalog No: #HW205

Package Size: #HW205-1 50ul #HW205-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2.00E+07
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7) detects endogenous levels of HIST1H3A
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Histone H3
Other Names	HIST1H3A
Accession No.	Swiss-Prot:P68431 Gene ID:8350
Uniprot	P68431
GeneID	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
Calculated MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl Histone H3 K9 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Stomach Tissue using Acetyl Histone H3 K9 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only