

Carbonic anhydrase 3

Catalog No: #AP81087

Package Size: #AP81087-1 50ug #AP81087-2 100ug #AP81087-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carbonate dehydratase III,Carbonic anhydrase III,CA-III,Car3
Accession No.	Uniprot ID:P16015
Uniprot	P16015
GeneID	12350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ca3

Note: This product is for in vitro research use only