

B-cell lymphoma/leukemia 11B

Catalog No: #AP81090



Package Size: #AP81090-1 50ug #AP81090-2 100ug #AP81090-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell lymphoma/leukemia 11B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BCL-11B,B-cell CLL/lymphoma 11B,COUP-TF-interacting protein 2,Radiation-induced tumor suppressor gene 1 protein,hRit1,CTIP2,RIT1
Accession No.	Uniprot ID:Q9C0K0
Uniprot	Q9C0K0
GeneID	64919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCL11B

Note: This product is for in vitro research use only