

A disintegrin and metalloproteinase with thrombospondin motifs 20

Catalog No: #AP81092

Orders: order@signalwayantibody.com

Package Size: #AP81092-1 50ug #AP81092-2 100ug #AP81092-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM-TS 20
Accession No.	Uniprot ID:P59511
Uniprot	P59511
GeneID	223838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adamts20

Note: This product is for in vitro research use only