

## Epsilon Tubulin Mouse Monoclonal Antibody(3G1)

Catalog No: #44172

Package Size: #44172-1 50ul #44172-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Epsilon Tubulin Mouse Monoclonal Antibody(3G1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3G1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Epsilon Tubulin Mouse Monoclonal Antibody(3G1) detects endogenous levels of Epsilon Tubulin
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Epsilon Tubulin
Other Names	Epsilon Tubulin
Accession No.	Swiss-Prot:Q9UJT0Gene ID:51175
Uniprot	Q9UJT0
GeneID	51175;
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only