

Uncharacterized protein C1orf106

Catalog No: #AP81100



Package Size: #AP81100-1 50ug #AP81100-2 100ug #AP81100-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Uncharacterized protein C1orf106
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q3KP66
Uniprot	Q3KP66
GeneID	55765;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C1orf106

Note: This product is for in vitro research use only